

## Production scale up of SDVB monolithic columns for the purification of mRNA molecules

T.Drožina<sup>1</sup>, P.Lapajne<sup>1</sup>, B.Košir<sup>1</sup>, A.Štrancar<sup>1</sup>, I.Bergoč<sup>1,\*</sup>

<sup>1</sup> BIA Separations d.o.o., a Sartorius Company, Mirce 21, 5270 Ajdovščina, Slovenia  
\* Corresponding author: ines.bergoc@biaseparations.com

### Introduction

Messenger RNA (mRNA) is becoming a major contributor in the fields of **gene therapy** and **vaccines**, including those developed in response to the COVID-19 pandemic. Convective Interaction Media<sup>®</sup> (CIM<sup>®</sup>) Styrene-divinylbenzene (SDVB) monolithic columns are promising for high resolution purification and separation of mRNA, enabling large-scale production of this molecule.

This study demonstrates the ability to prepare homogeneous SDVB monoliths with desired chromatographic properties and economical analytics over the whole size range.

### 1. Experimental approach

SDVB monoliths were prepared by free-radical thermal polymerization of styrene and divinylbenzene under isothermal reaction conditions. The desired pore-size distribution was achieved using conventional porogenic solvents. The prepared SDVB columns exhibited high degree of crosslinking with desired morphology. A well-controlled radical polymerization process was transferred from the laboratory to industrial scale, i.e., from 1 mL columns up to 8000 mL CIMmultus<sup>™</sup> columns. Furthermore, the production of SDVB monoliths was scaled down to 0.1 mL CIMac<sup>™</sup> columns, demonstrating our ability to control the polymerization process over the whole size range.

### 2. Gating Strategy

Another challenge in this study was to develop inline analytics that could confirm monolith homogeneity and desired chromatographic properties for columns from 0.1 to 8000 mL. For this purpose, two HPLC chromatographic methods were used. The cGMP compliant is an isocratic one, separating unretained uracil and retained small molecule, while the bioapplication-based employs size separation of two messenger RNA molecules in an acetonitrile gradient.




Figure 1: Diagram of gating strategy of production scale up and correlation between cGMP compliant and bioapplication-based method

### 3. Results – Method for Test Chromatogram

An isocratic liquid chromatographic method with UV detection at 270 nm was used for standard test chromatogram (TK) to verify packaging conditions and column performance. Mobile phase consisted of a mixture of ethanol and salt. The sample used consisted of an unretained uracil and a retained small molecule. The loop volume and flow rate were set according to column size. To eliminate the influence of different HPLC systems two chromatographic parameters were considered, i.e., retention volume of HPLC system with column ( $V_r$ ) and normalized retention factor of sample ( $k'$ ).




Figure 2: Typical test chromatogram of 80 mL column with retained small molecule and unretained uracil for monitoring chromatographic properties and repeatability of produced SDVB columns.

### 4. Results – Method for mRNA molecules

The gradient separation method with acetonitrile was used to test performance of 0.2 mL columns, which were cut from the large monolith. The normalized retention factors of two mRNA molecules ( $K_{1000\text{ nt}}$  and  $K_{4000\text{ nt}}$ ) were calculated.




Figure 3: Representative mRNA chromatogram on 0,2 mL specimen of two mRNA molecules

### 5. Results – Correlation between TK and mRNA





Figure 4: Correlation between normalized retention factor of TK and mRNA chromatographic methods.

$$k' = \frac{t_{\text{retained comp.}} - t_{\text{unretained comp.}}}{t_{\text{unretained comp.}} - t_{\text{SST}}}$$

$$K_{1000\text{ nt}} = \frac{t_{1000\text{ nt}} - t_{\text{unretained comp.}}}{t_{\text{unretained comp.}}}$$

$$K_{4000\text{ nt}} = \frac{t_{4000\text{ nt}} - t_{\text{unretained comp.}}}{t_{\text{unretained comp.}}}$$



### 6. Conclusion

- Correlation between cGMP and bioapplication-based methods
- Scale up from 0,1 mL to 8000 mL columns
- highly reproducible