

July 19, 2023

Keywords or phrases:

Adeno-associated virus, capture step, chromatography, monolith, purification, optimization

Fast and Semi-automatic Optimization of Capture Step for rAAV Purification Using CIM® SO3 96-Well Monolithic Plate

Kaja Bažec, Mirjam Krašna, Andrej Mihevc, Maja Leskovec, Aleš Štrancar, and Mojca Tajnik Sbaizero*
Sartorius BIA Separations d.o.o., Mirce 21, 5270 Ajdovščina, Slovenia

* Correspondence
Email: mojca.tajniksbaizero@sartorius.com

Abstract

Adeno-associated viruses (rAAVs) are commonly used vehicles for gene therapy applications, leading to increased demand for scalable purification tools. However, since rAAVs are a relatively new and highly variable modality, established platform processes do not yet exist. This can create long process development timelines.

Here, we show that CIM® SO3 monolithic 96-well plates can be used to screen cation exchange capture chromatography parameters and that results are scalable to preparative scale on the CIMmultus® line, which is suitable for industrial-scale applications.

Introduction

A robust capture step is a key parameter to increase yield during rAAV production. During monolithic chromatography capture, the rAAV is selectively bound to a strong cation exchanger (SO3) in acidic pH conditions.^{1,2} In this step, the majority of host cell contaminants are removed, and rAAV particles are strongly bound to the matrix. The productivity of the step is strongly influenced by the dynamic binding capacity (DBC), which depends on the sample preparation and the availability of the binding sites at the chromatographic matrix.

As the demand for rAAV quantity and quality in clinical manufacturing continues to rise, effective solutions for industrial-scale purification of large volumes of upstream material are required. rAAV capture is the workhorse of the chromatography workflow and presents a valuable optimization opportunity through which the full potential of SO3 monolithic columns can be realized.

To reduce the development time using multi-factor screening and increase comparability during downstream process development, CIM® SO3 monolithic 96-well plates were introduced to the Sartorius BIA Separations portfolio. Here, we show that results on this small scale can be applied to CIMmultus® columns, which is scalable to large volumes.

Methods and Results

Fine-Tuning Mobile Phase Composition Using CIM® SO3 96-Well Plate Screening


CIM® SO3 monolithic 2 µm 96-well plates were used to optimize mobile phase screening for rAAV downstream process development. Using a vacuum pump system, each well was prefilled with 0.05 mL of 2 µm SO3 monolith to ensure uniform flow and high flow rates across the plate.

To optimize mobile phase A, we screened buffers of different pH and NaCl concentrations. Each buffer combination was also tested in the presence of Poloxamer 188 to determine whether this additive influences the binding of the product to the SO3 monolith. Each combination was tested in triplicate (Figure 1). The elution buffer (mobile phase B) contained 2 M NaCl and the same pH and Poloxamer 188 concentration as mobile phase A for each condition.

Before loading on the SO3 96-well monolithic plate, clarified rAAV9 lysate obtained from Sf9 cells was processed by tangential flow filtration (TFF) coupled with nuclease treatment – salt-tolerant nuclease, as previously described.³ The TFF retentate was acidified to the selected pH (pH 3.5, 4.0, 4.5, or 5.0) and filtered using Minisart® PES 0.45 µm filters.

After removing the storage solution and washing with diH₂O, the CIM® SO3 96-well plate was conditioned with mobile phase A for 10 CV (CV; column volume, 0.05 mL), mobile phase B for 10 CV, and mobile phase A once more for 15 CV. The wells were loaded with acidified filtered TFF retentate (1 mL), and the flowthrough fraction was collected into a collecting plate. Next, the column was washed with 10 CV of tested mobile phase A, and the wash fraction was collected into the collecting plate. Elution was performed using 4 CV of mobile phase B and also collected. CIP was performed with 10 CV of CIP solution (1 M NaOH 2 M NaCl). After this, the plate was regenerated and stored per the manufacturer’s instructions.

Figure 1: Design of SO3 Monolithic 96-Well Plate Testing




Automating Analytics for CIM® SO3 96-Well Elution Plate Samples

One of the main advantages of using 96-well monolithic plates is the standard 96-well design compatible with all appropriate 96-well analytical methods. This saves a significant amount of time, is user-friendly, and reduces the number of steps needed to gain information about the viral titer and impurities for each of the tested conditions during process development.

The wash and elution fractions for mobile phases A and B were analyzed using a fluorescence plate reader with Tryptophan native fluorescence settings (ex. 280 nm, em. 348 nm). This allowed us to follow where concentrated AAV particles eluted under each condition (Figure 2). The results show that higher pH is negatively correlated with increased NaCl concentration in the wash buffer (mobile phase A). For pH 3.5, all tested NaCl conditions up to 750 mM showed binding of AAV to the SO3 monolith (high fluorescence detection values in elution fractions), whereas binding at higher pH resulted in tolerating lower NaCl conditions.


Figure 2: Fluorescence Readings (A) Wash (Mobile Phase A) and (B) Elution (Mobile Phase B) Fractions



To test the efficiency of rAAV9 binding and elution, capsid-specific ELISA was performed for the conditions that were close to maximum NaCl for each tested pH, as well as predicted negative binding conditions (low fluorescence signal in elution fraction). Viral titer (vp/mL) was the highest for mobile phase A, pH 3.5, and 500 mM NaCl in the presence of Poloxamer 188 (Figure 3A). To test the impurity profile, we performed SDS-PAGE analytics in the same elution fraction as ELISA was performed. Qualitatively, there is an observed trend towards increased impurities in the load (acidified TFF retentate) sample. More low-molecular-weight impurities were also observed in elution fractions from high pH samples, compared to lower pH elution (Figure B).

The combination of parameters (different pH, sodium chloride, and poloxamer concentration) that gave us the highest virus titer (vp/mL) value and lowest impurity content in elution fractions was selected as the most optimal mobile phase A (pH 3.5 and 500 mM NaCl in the presence of Poloxamer 188).


Figure 3: (A) rAAV9 Viral Titers (vp/mL; (ELISA) and (B) Impurity Fractions (SDS-PAGE) in Elution Fractions



Scale-Up to CIMmultus® Line

We then assessed the scalability of these results at preparative scale. We tested selected chromatographic conditions from the CIM® SO3 96-well monolithic plate on a 1 mL CIMmultus® SO3 column. We applied a linear NaCl gradient to observe the rAAV elution profile (Figure 4A). The results show that the capture step was robust under the selected conditions; high purity was observed in the elution AAV fraction (Figure 4B). With AAV9-capsid-specific ELISA, we detected 87.7% rAAV9 capsids in the E2 (main peak) fraction compared to the SO3 load material.

Figure 4: (A) Preparative Chromatogram Using CIMmultus® SO3 1 mL Column (Conditions Selected From CIM® SO3 96-Well Plate Screening); (B) SDS-PAGE Analytics for Tested Fractions.



Conclusion

CIM® 96-well plates are a useful tool for fast and automated screening of parameters during process development. A low sample amount is required for experiments, and several parameters can be studied in parallel to achieve high comparability. As such, CIM® SO3 96-well plates enable fine-tuning of an rAAV capture step, which can be transferred and adapted to the CIMmultus® SO3 product line.

References

1. HEK293-Derived Adeno Associated Virus (AAV) Purification: Comparison of Small Scale Laboratory Production Towards Industrial Format Using Monoliths. Poster, BIA, <https://www.biaseparations.com/en/library/posters/1002/hek293-derived-adeno-associated-virus-aav-purification-comparison-of-small-scale-laboratory-production-towards-industrial-format-using-monoliths>
2. An Industrial Purification Platform for AAV. Poster BIA, <https://www.biaseparations.com/en/library/posters/988/an-industrial-purification-platform-for-aav>
3. Gagnon P., Leskovec M., Goričar B., Štrancar A. Streamlining Industrial Purification of Adeno-Associated Virus. Bioprocess International, December 2020.

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen
Phone +49 551 308 0

USA

Sartorius Stedim North America Inc.
565 Johnson Avenue
Bohemia, NY 11716
Toll-Free +1 800 368 7178

Slovenia

Sartorius BIA Separations
Mirce 21
5270 Ajdovscina
Toll-Free +386 59 699 500

 **For more information, visit**
www.sartorius.com