



A045

Purification of influenza A virus from embryonated chicken eggs using CIMac™ Adeno-0.1 Analytical Column

Influenza vaccines are still predominantly produced in embryonated chicken eggs and the purification processes barely have changed during the years. There is a growing need for fast, efficient and economical vaccine production.

So far, monolithic supports have been used successfully in virus purification and concentration, as well as in the purification of virus-like particles (VLP) propagated in cell cultures.

Therefore, our aim was to prove the applicability of monoliths in purification of influenza virus A propagated in embryonated chicken eggs.



SAMPLE PREPARATION

Influenza A virus (equine H3N8 strain) was inoculated into the allantoic cavity of 10 day old embryonated chicken eggs.



RESULTS


> Figure 1

Step gradient elution of influenza virus harvest

Column:	CIMac™ Adeno-0.1 Analytical Column (Pores 2 µm); Catalog number: 110.8502-2
Column chemistry:	Strong anion exchanger; quaternary amine (QA)
Instrumentation:	HPLC Knauer system
Mobile phases:	Buffer A (MPA): 50 mM HEPES; pH 7.5 Buffer B (MPB): 50 mM HEPES containing 2 M NaCl; pH 7.5
Loading material preparation:	The virus harvest was diluted 10 times in equilibrating mobile phase (MPA) and filtered through 0.45 µm CA membrane filters.
Column loading volume:	1.32 mL
Flow rate:	1 mL/min
Method:	Step gradient elution: 0.2 M – 0.5 M – 0.8 M – 2M NaCl in MPA
UV detection:	UV at 280 nm
Virus quantification:	HA assay
Protein quantification:	Bradford Ultra assay
DNA quantification:	PicoGreen assay

> **Figure 2**

Step gradient elution chromatogram of influenza virus on CIMac™ Adeno-0.1 Analytical Column



> **Table 1**

Recovery of influenza virus, proteins and DNA during the purification process

Fraction	Virus's HA titer (% of load)	Protein (% of load)	DNA (% of load)
Load	100	100	100
FT + W	0	15	5
E1	0	68	23
E2	88	13	10
E3	6	0	42
E4	0	0	4
Total recovery	94	96	84

CONCLUSIONS:

An efficient purification process for influenza A virus obtained from allantoic fluid was developed on CIMac™ Adeno-0.1 Analytical Column (quaternary amine). High virus yield (88%) and efficient depletion of proteins (87 % total proteins) and DNA (90% total DNA) was achieved.

CIM® technology enables fast and efficient one step purification of influenza A virus produced in embryonated chicken eggs.



For any additional information please contact us:

Tel.: +386 5 9699 500

sales@biaseparations.com

www.biaseparation.com

Information and specifications contained here are, to the best of our knowledge, accurate and represented in good faith. They are intended to help you start working with this new separation technology and are subject to change without notice. BIA Separations shall not be liable for errors contained herein or for incidental or consequential damages in connection with the performance of use CIM. For more information on our products, visit our home page at: <http://www.biaseparations.com> or contact your local distributor. We reserve the right to alter the specification detail etc. without prior notice or liability.