

# SARTORIUS

Simplifying Progress



Orthogonal PAT analytics - essential to enable robust production of viral vectors and LNPs

Aleš Štrancar

Sartorius BIA Separations, Ajdovščina, Slovenia

# Sartorius BIA Separations product and solutions portfolio

## BIA's Proprietary Product and Services Offering

### CIMmultus Monolith production-scale columns



### PATfix systems & CIMac Monolith process control columns



### Process and method development services



**10 - 100x**  
higher capacity

**up 3x**  
higher recovery

Higher integrity information  
at lower cost and  
better accuracy

Achieves rapid,  
high-resolution  
separations in minutes

The leading expert on large  
biomolecules and viral  
particles

Tailored services for each  
clinical phase up to and  
including production

**up to 3x**  
lower drug  
manufacturing cost

Column sizes  
**40L+**

Biocompatible

Fast and reproducible LC  
monitoring of  
large biomolecules

Drives long-term,  
embedded customer  
relationships

50+ Process and Method  
development experts

**Serving Gene Therapy for > 25 years**

# Experts in fast bioproduct manufacturing process development – no royalties charged

>80 pDNA, mRNA, virus DSP  
cGMP processes tech transferred  
to CMOs, sponsors, including  
Corona.

Product impurities are one of the  
key reasons for treatment side  
effects. High purity is therefore  
mandatory for product safety.

- **pDNA, purity is THE key for better transcription and better transfection for purer mRNA**
- Minicircle DNA, ministring (shorten the pDNA), ssDNA
- **ssRNA and dsRNA, platform process from E.coli to mRNA and LNP**
- Adeno virus, more than 20 years experience
- **AAV (all serotypes, > 20 tested)**
- Life influenza virus (all serotypes)
- Vaccinia/MVA
- Exosome
- Bacteriophages (outstanding endotoxin removal)
- VLPs and inactivated vaccines
- **LNPs (single, multiple payloads, “decorated”)**
- IVIG
- IgM and more

# Testimonials: need for speed

*“We are especially grateful that BIA Separations shared, and operated, with the same sense of urgency we did to help bring gene therapy to the SMA community. BIA’s experience with AAV purification and its chromatographic technology were important contributions and we look forward to our continued work together.”*

*Andy Stober, Senior Vice President of Technical Operations for AveXis/Novartis*

*“We use BIA Separations Monolith columns for the purification of our mRNA drug substance. The Monolith columns come in multiple sizes to meet our needs from small scale product development work to large scale cGMP manufacturing runs. We have found the Monolith columns to provide high throughput and high purity while being very robust and reliable. In addition, BIA Separations has top notch customer service that never fails to impress.”*

*Greg Kubczak, Director of Technical Services and Manufacturing at Arcturus Therapeutics*

*“We are very pleased to be collaborating with BIA Separations / Sartorius in developing and tech-transferring optimized mRNA vaccine manufacturing processes. When time is of the essence, it is crucial to be able to rely on key partners whilst maintaining the highest level of quality in every aspect of our work. We are also very proud that our team was able to produce two kinds of genetic vaccines (DNA and mRNA) and we look forward to expanding our partnership with the BIA Separations’ team on future projects.”*

*Hong Thai Pham, CEO at BioNet*



<https://www.zolgensma.com>

# KEY MESSAGES

## **Do NOT trust analytics**

Orthogonal analytical methods are mandatory to prevent wrong conclusions due to the analytical method misleading, and should be in place before you start process development

PATfix systems were developed to address these requirements

**Process is as good as your analytics are**  
**Product is as good as your analytics are**

**The more defined the process, the easier it's scale-up would be.**

# PATfix, all in one system, provides orthogonality to molecular biology methods



- Analyses of various biomolecules; pDNA, mRNA, dsDNA, protein...



- Used for process monitoring or release of product and control



- Rapid chromatographic analytical methods (~3-15 min)



- HPLC analytics does not require expensive reagents



- Broad diversity of column chemistries; PrimaS, Oligo-dT, pDNA, SDVB...



- Various analytical techniques; pH, gradients, RP, SEC, HIC

## Qualified/validated methods include:

- AAV particle quantification
- AAV empty/partial/full determination
- Aggregate quantification
- Residual hc DNA and hc proteins
- Empty/full in the bioreactor and more

## NO need to...



- Search for the right LC hardware



- **Selecting suitable column**



- Develop an internal standard



- **Develop/optimize the analytical method**




- Validate the method



# Fast process development enabled by multidetector PATfix

**PATfix with multiple detectors allows for sample characterisation in an hour.**

**Novel, very reliable and affordable MLS detector, including the prep version.**



PATfix and CIMac →  
Process understanding

**Cation exchange** does not discriminate empty from full capsids but it still provides fast characterization of total AAV and contaminant content. UV wavelength ratios provide a hint about relative DNA and protein distribution but **fluorescence enables direct quantitative comparison**.



From *E.coli* to LNP with one dimension PATfix system

Simplifying Progress




**SARTORIUS**


# From *E. coli* to mRNA

## Sartorius mRNA Drug Substance Production Workflow

Analytical workflow (PATfix pDNA; CIMac pDNA)



E.Coli Fermentation



Alkalizator  
CaCl<sub>2</sub> treatment



UF | DF



CIMmultus  
DEAE



Linearization



CIMmultus  
C4 HLD



UF | DF




Linear  
pDNA



Linear  
pDNA



IVT in AMBR250  
reactor




CIMmultus Toolbox  
Oligo dT, PrimaS, C4 HLD, SDVB



UF | DF



mRNA



LNP

Analytical workflow (PATfix mRNA; CIMac PrimaS, CIMac Oligo dT, CIMac SDVB)

# PATfix, all in one system, provides orthogonality to molecular biology methods


-  ▪ Analyses of various biomolecules; pDNA, mRNA, dsDNA, protein...
-  ▪ Used for process monitoring or release of product and control
-  ▪ Rapid chromatographic analytical methods (~3-15 min)
-  ▪ HPLC analytics does not require expensive reagents
-  ▪ Broad diversity of column chemistries; PrimaS, Oligo-dT, pDNA, SDVB...
-  ▪ Various analytical techniques; pH, gradients, RP, SEC, HIC

NO need to...

-  ▪ Search for the right LC hardware
-  ▪ Selecting suitable column
-  ▪ Develop an internal standard
-  ▪ Develop/optimize the analytical method
-  ▪ Validate the method

## Qualified/validated methods include:

- pDNA purity and quantification
- mRNA IVT monitoring
- mRNA quantification and purity
- **LNP encapsulation, purity, size, heterogeneity, lapidated RNA amount,**
- **LNP – protein conjugation reaction monitoring**
- **and more in development**





Two dimensions PATfix LNP  
Switcher for LNP process  
development and production


Simplifying Progress




**SARTORIUS**

# Encapsulation efficiency analysis using PATfix LNP Switcher (two dimensional LC)




Direct injection of LNP formulations on LC Shearless monolithic column; no sample pre-treatment, dilution with loading buffer only



# Encapsulation efficiency analysis using PATfix LNP Switcher – chromatogram



# PATfix analysis orthogonal to LNP size & LNP size distribution measurements



- OH column separates particles into populations by size
- LC LNP Switcher enables determination of the size of lipid nanoparticles by utilising the MALS detector
- Radius was determined by modelling 9 scattering angle responses using the Berry model and a 40 nm PS standard
- LC method in agreement with Nano Tracking Analysis (NTA)

# Encapsulation efficiency using PATfix LNP Switcher analysis – orthogonal to RiboGreen

## Quant-it™ RiboGreen method

- RiboGreen dye is used for detection and quantification of nucleic acids
- Triton X-100 is applied for LNPs opening
- There is no significant difference in the results obtained with these two techniques


Comparison of encapsulation efficiency results obtained by LNP Switcher and Quant-it™ RiboGreen assay:

Sample	% encapsulation	
	LC LNP Switcher	Quant-it™ RiboGreen
LNP mRNA	72	71




PATfix method can be used as orthogonal to RiboGreen method

# PATfix LNP Switcher analysis of co-encapsulated cargo – LNP homogeneity analysis




- LC switcher separates and quantifies co-encapsulated species, such as Cas9/sgRNA for CRISPR applications.
- Separate value of encapsulation efficiency for each mRNA can be assessed.
- **LNP can have very high encapsulation efficiency, but can be very heterogeneous.**

# PATfix LNP Switcher allows for monitoring of lipidated mRNA impurities



- Enabled monitoring of mRNA-lipid covalent adducts formed from ionizable lipid impurities.
- Important for lipid screening and stability studies.
- Highly dependent on lipid choice, storage conditions, and raw material purity.
- LNP Switcher, unlike CE, does not need any sample pre-treatment!

# PATfix LNP Switcher for CRISPR formulations



Just one sample injection allows for simultaneous:

- Separation and quantification of multiple-cargo formulations (e.g. CRISPR, CAR-T, combination vaccines)
- Accurate quantification of different nucleic acid species
- Encapsulation efficiency determination
- Heterogeneity assessment
- RNA-lipid adduct quantification
- RNA fragmentation assessment



Targeted LNP (tLNP) – conjugation  
monitoring

Simplifying Progress




**SARTORIUS**

# Conjugation reaction monitoring



- On-particle conjugation reaction (incubation with protein for 2h)
- Two chromatographic methods were used to monitor conjugation reaction process:  
HIC (CIMac OH column - free protein monitoring) and AEX (CIMac QA column (targeted vs untargeted LNPs separation))
- Reaction progress tracked at different times after incubation starts:
  - t = 0 min
  - t = 20 min
  - t = 40 min
  - t = 60 min
  - t = 120 min
  - t = 1 day

# Conjugation reaction monitoring - PATfix CIMmultus OH




- t = 0 (red): untargeted particles, no proteins in flowthrough
- With progressing reaction, (t)LNP peak in gradient is shifting to the left, area is getting larger
- Fluorescence peak area is getting bigger as reaction progresses (more protein on the surface), but MALS peak stays roughly the same


Zoom-in



# Orthogonal conjugation reaction monitoring - PATfix CIMmultus QA HR



- t = 0 (red): untargeted particles, no FT
- Peak in gradient starts to shift to the right, fluorescence signal is getting bigger (more protein on particles), but MALS signal stays roughly the same




# Conjugation reaction monitoring using PATfix LNP switcher

➤ LNP switcher used to simultaneously monitor:

- the conjugation reaction between the protein and the LNP particle
- nucleic acid content of LNPs
- LNP encapsulation efficiency, concentration, lipidated RNA,...

➤ Monitoring the conjugation reaction is the key to optimize formulations performance and control the manufacturing process





Shearless LNP manufacturing

# Shearless chromatography to replace TFF degrading LNP products

**Challenge**




Avoid scale-up issues and shear forces.  
Concentrate LNPs and remove free mRNA.  
Retain uniform and functional particles.  
Remove EtOH





**Solution**



Low shear swift preparatory hydrophobic interaction chromatography on monolithic columns, CIMmultus OH



# Shearless chromatography improves the potency / recovery



Purifying LNPs by chromatography using CIMmultus OH separates particles by size and can be used to remove subpopulations with lower potency.

# Conclusions

- Analytical methods might be misleading. Orthogonal methods are needed for correct process understanding and its control.
- LNP PATfix Swither can be used as orthogonal to RiboGreen method and to DLS/NTA/Videodrop methods, and in addition allows for:
  - LNP encapsulation,
  - LNP purity,
  - LNP heterogeneity,
  - Iapitated RNA amount,
  - LNP – protein conjugation reaction monitoring, and
  - LNP corona studies (coming).
- Purification using CIM OH monolithic columns gives recoveries >90%, even of more complex LNPs, and produces particles of higher uniformity and potency.
- Orthogonal in-process analytics, THE key pillar for the CMC.

# KEY MESSAGES

Do NOT trust analytics

Orthogonal analytical methods are mandatory to prevent wrong conclusions due to the analytical method misleading, and should be in place before you start process development

PATfix systems were developed to address these requirements

Process is as good as your analytics are  
Product is as good as your analytics are

The more defined the process, the easier it's scale-up would be.

# Thank you for your attention!

Visit us at the Sartorius BIA Separations  
booth at exhibition hall

[ales.strancar@biaseparations.com](mailto:ales.strancar@biaseparations.com)

**SARTORIUS**