






Fichtenbaum, Andreas 1; Pungor, András 2; Nagy, Zoltán 2; Černigoj, Urh 3; Vidič, Jana 3; Barut, Miloš 3; Mitulović, Goran 1  
 1Medical University of Vienna; 2University of Debrecen, Debrecen, Hungary; 3BIA Separations d.o.o., Ajdovščina, Slovenia

## INTRODUCTION

Enrichment of phosphopeptides prior to LC-MS analysis is a crucial sample preparation step because of their low stoichiometry in biological sample, longer retention on reversed phase columns, and lower ionization efficiency compared to non-phosphorylated peptides [1]. The use of metal oxides, most prominently of TiO<sub>2</sub> enabled efficient and relatively simple phosphopeptide-enrichment. In this study a new monolithic column from BIA Separations containing immobilized TiO<sub>2</sub>-nanoparticles was tested for its ability to enrich phosphopeptides. The TiO<sub>2</sub>-column was also tested for possible carryover originating from biological samples. In conclusion, tested monolithic TiO<sub>2</sub> columns show significant binding ability for phosphopeptides and are considered as suitable for phosphopeptide enrichment.

## MATERIAL AND METHODS



## RESULTS - Phosphopeptides retained on TiO<sub>2</sub> column versus Blank column




Fig.4: Comparison of identified phosphopeptides of tryptical digested bovine α-casein in the elution fraction "E1" collected from phosphopeptide enrichment on a monolithic column with and without immobilized nano-TiO<sub>2</sub> particles.

## RESULTS - Carryover of the monolithic nano-TiO<sub>2</sub> column




Fig.7: Testing carryover of the nano-TiO<sub>2</sub> monolithic column after phosphopeptide enrichment with α-casein (brown) and a blank sample (blue).

## RESULTS - Phosphopeptides enriched on the monolithic nano-TiO<sub>2</sub> column




Fig.5: Number of phospho- and non-phosphorylated peptides from tryptically digested bovine α-casein. Fractions (FT,W1,W2,E1,E2) were collected upon phosphopeptide enrichment on a monolithic column containing immobilized nano-TiO<sub>2</sub> particles. α-casein digest (Digest) directly injected to LC-MS served as control sample.

Fig.6: Phosphopeptides from α-casein identified in the elution fraction "E1" after phosphopeptide enrichment on the monolithic nano-TiO<sub>2</sub> column.

## DISCUSSION

A new monolithic column containing superficially immobilized TiO<sub>2</sub> nano-particles was tested for its ability to enrich phosphopeptides from a tryptical digest of bovine α-casein. Binding ability of phosphopeptides to the particles was confirmed compared to a blank column (Fig.4).

As expected the majority of non-phosphorylated peptides didn't bind to the nano-TiO<sub>2</sub> support. These were detected in the flow-through (FT) fraction. Elution fraction E1 contained most unique phosphorylated peptides (Fig.5).

Results of this evaluation study indicate the suitability of the tested nano-TiO<sub>2</sub> coated monolithic column for enrichment of phosphopeptides.

Carryover on the monolithic nano-TiO<sub>2</sub> column was not observed except in elution fraction "E1" the phosphopeptide VPQLEIVPNSAEER was identified (Fig.7). Since this phosphopeptide was found constantly in all fractions (FT,W1,W2,E1,E2) it indicates a strong binding affinity to the column. Therefore, further protocol optimization in terms of column cleaning is suggested.

## CONCLUSION

Mesoporous nano-TiO<sub>2</sub> particles immobilized on a monolithic support may provide higher capacity than common TiO<sub>2</sub> material and lower backpressure in HPLC which makes it capable also being applied to a nano-sized approach. Evaluating the binding ability and carryover of the column was a first step in this study but protocol optimization has to be evaluated especially on tryptical digests derived from a complex biological sample.

Moreover, since a study of Bodenmiller et al [3], that compares several phosphopeptide enrichment methods revealed only 30% of overlapping phosphopeptides, another outlook to further improve MS identification rates on phosphorylations is the combination of the nano-TiO<sub>2</sub> support with a second enrichment technique such as e.g. Electrostatic Repulsion Chromatography.

## REFERENCES

1. L. Beltran et al, Amino acids 43, 1009 (Sep, 2012).
2. M. Mazanek et al., Journal of chromatography. B, Analytical technologies in the biomedical and life sciences 878, 515 (Feb 15, 2010).
3. B. Bodenmiller et al, Nature methods 4, 231 (Mar, 2007).